

Insert racks for Waterbaths

Insert rack for test tubes and butyrometers

Polypropylene (PP) white, acid and alkaline-proof, heat-resistant up to 100°C Levels connected by sleeve screws (Height dimension with handle)



Order Details:

Quantity	Drilling-D/mm	LxWxH/mm	Order No.
12	26	220x80x220	ESGRB 12
32	19	224x125x147	ESGRB 32
56	17	380x120x220	ESGRB 56

Insert rack for Mojonier Pipes

Polypropylene (PP) white, acid and alkaline-proof, heat-resistant up to 100°C Levels connected by sleeve screws (Height dimension with handle)



Order Details:

Quantity	Drilling-D/mm	LxWxH/mm	Order No.
8	39	380x100x220	ESGM

Insert rack for Mikli and baby bottles

Polypropylene (PP) white, acid and alkaline-proof, heat-resistant up to 100°C Levels connected by sleeve screws (Height dimension with handle)



Order Details:

Quantity	Drilling-D/mm	LxWxH/mm	Order No.
----------	---------------	----------	-----------

12 (2×6)

54

380x120x220 ESMB
