


Racks & Supports

SmartLab Racks with Numbering

The Aerne Analytic SmartLab Series includes test tube racks for 6, 12, 24, 36 and 48 samples for a diameter of each 18mm. Compared to conventional test tube racks the sample spaces of the SmartLab Series is graved with numberings to facilitate the daily work in laboratories. There is no more lettering of the single samples necessary. The racks are made of PVC and therefore are perfect for the use in laboratories because of their stability and possibility to clean it sterile.




Order details:

Part	Order-No.
for 12 (2×6) glasses	SLS 6
for 12 (2×12) glasses	SLS 12
for 24 (3×12) glasses	SLS 24
for 36 (4×12) glasses	SLS 36
for 48 (4×12) glasses	SLS 48

Racks for butyrometers

Polypropylene (PP) white, acid and alkaline-proof, heat-resistant up to 100°C, Drillings 26 mm diameter



Order details:

Part	Order-No.
for 12 (2×6) glasses; Dimensions (mm): L=160 x W=60 x H= 87	RGG 12
for 24 (2×12) glasses; Dimensions (mm): L=315 x W=60 x H= 87	RGG 24

for 36 (3×12) glasses; Dimensions (mm): L=315 x W=90 x H= 87 RGG 36

for 48 (4×12) glasses; Dimensions (mm): L=315 x W=120 x H= 87 RGG 48

Racks for butyrometers

Polypropylene (PP) white, acid and alkaline-proof, heat-resistant up to 100°C, Drillings 26 mm diameter




Order details:

Part	Order-No.
for 12 (2×6) Butyrometer; Dimensions (mm): L=225 x W=80 x H= 155	BMG 12
for 24 (2×12) Butyrometer; Dimensions (mm): L=450 x W=60 x H= 155	BMG 24
for 36 (3×18) Butyrometer; Dimensions (mm): L=630 x W=90 x H= 155	BMG 36

Rack for centrifuge tubes (ZG)

Polypropylene (PP) white, acid and alkaline-proof, heat-resistant up to 100°C



Order details:

Part	Order-No.
Number centrifuge tubes: 24(2×12) ;Drilling: 18 mm; L:315xW:60xH:65	ZGG 2/12/18
Number centrifuge tubes: 12(2×6); Drilling: 26 mm; L:220xW:80xH:65	ZGG 2/6/26

Number centrifuge tubes: 6 (1×6) ;Drilling: 36 mm;
L:315xW:60xH:65

ZGG 1/6/36

Number centrifuge tubes: 6 (1×6) ;Drilling: 42 mm;
L:315xW:60xH:65

ZGG 1/6/42

Imhoff cone Racks

PVC white, rack rods 500mm long, height-adjustable, Drillings 98 mm diameter



Order details:

Part

Order-No.

Number of cones: 2; Size of Base plate: 310x150mm (Höhe 500mm)

SIG 2

Number of cones: 4; Size of Base plate:650x150mm (Höhe 500mm)

SIG 4

Number of cones: 6; Size of Base plate: 900x150mm (Höhe 500mm)

SIG 6

Racks for Funnels

PVC white, rack rods 500mm long, height-adjustable, Drillings 60 mm diameter



Order details:

Part

Order-No.

Number of Funnels:2; Size of Base plate: 230x100mm (Höhe 500mm)

TGS 2


Number Trichter:4; Size of Base plate:500x100mm (Höhe 500mm)

TGS 4

Number Trichter:6; Size of Base plate:770x100mm (Höhe 500mm)	TGS 6
Number Trichter:10; Size of Base plate:1040x100mm (Höhe 500mm)	TGS 10

Rack for Araeometers

PVC white, with drip tray, for convenient storage of Araeometers up to 30 mm diameter




Order details:

Part	Order-No.
Rack width 140mm (Höhe 200mm) for 4 Araeometers	AMG 4
Rack width 206mm (Höhe 200mm) for 6 Araeometers	AMG 6

Racks for weighing pipettes

PVC white, for bent weighing pipettes



Order details:

Part	Order-No.
Rack for 2 Wägepipetten; Dimensions 65x35mm (Höhe 35mm)	WPG 2
Rack for 4 Wägepipetten; Dimensions 65x55 mm (Höhe 35mm)	WPG 4

Pipette Carousel



Order details:

Part

Pipette Carousel

Order-No.

PPK
