

Shaking Shells, Pour-out Caps, Shaking hood

Pour-out Caps for Butyrometers

PVC white

The pour-out cap is placed over the uncorked butyrometers and the rack which are held in place by the cap when poured or rinsed.




Order Details:

Part	Order No.
12 Test portions	AHB 12
24 Test portions	AHB 24
36 Test portions	AHB 36

Shaking hood for Butyrometers

PVC white

The shaking hood is placed over supports with sealed butyrometers. The side facing your body is protected by a dropping wall.



Order Details:

Part	Order No.
12 Test portions	SHB 12
24 Test portions	SHB 24
36 Test portions	SHB 36

Shaking shells for Butyrometers

Polypropylene (PP) white

The individual butyrometers are placed in the shell. Up to four butyrometers can be shaken, when the cork is held in place simultaneously.



Order Details:

Part	Order No.
1 Butyrometers	SHÜB 1
2 Butyrometers	SHÜB 2
3 Butyrometers	SHÜB 3
4 Butyrometers	SHÜB 4
