

Incubator Heatsensor



Application:

The procedure of inhibitor tests for milk analysis is one of the essential process steps for the quality assurance in the dairying. Test laboratories and farmers use the efficient quick test for their analysis to evaluate the today's high safety standards in this industry. The operating mode of this heating instrument embraces the constant tempering and incubation of the milk test. The Aerne Analytic incubators are equipped with aluminium heating units and an analogue/ digital temperature control.

Technical Data:

Controller:	Microprocessor
Number Test:	2 rows with 8 + 1 row with 6
Temperature:	40C +/- 0,30C
Timer:	3 + 3 Minutes with Signal Sound
Power Supply:	230 Volt AC, 110 Volt AC, 24 Volt DC, 12 VDC

Order Details:	
Digital:	Order-No.
230 Volt	2R-230D
12 Volt	2R-12D
24 Volt	2R-24D
110 Volt	2R-110D

Wärmeverteilung Aluminiumkörper

