

Cooling-Incubator

**Application:**

The Cooling-Incubator can be used for tempering of Samples at constant Temperature of 20°C (+/- 2°C). It is provided to be used in a dry labor enviroment.

Technical Data:

Housing:	PVC white
Controlling:	Foil Keyboard - JUMO
Adjustable Range:	15°C - 35°C +/- 1°C
Control:	JUMO - Quantrol LC100
Safety:	Shutdown of the Heating-Power at Overheating 77°C < 70 dB
Heater:	Heating - Resistor
	Material: CeramiC Power: 2 x 30 Watt
Kühlung:	Peltier Element: - Power: ca. 50 Watt
Power Supply:	230 VAC 50/60 Hz

Order Details:**Order-No.**

TT28B-230D