

Incubator CH-ATK for Chr. Hansen-Test



Anwendung:

The procedure of inhibitor tests for milk analysis is one of the essential process steps for the quality assurance in the dairying. Test laboratories and farmers use the efficient quick test for their analysis to evaluate the today's high safety standards in this industry. The operating mode of this heating instrument embraces the constant tempering and incubation of the milk test. The Aerne Analytic incubators are equipped with aluminium heating units and an analogue/ digital temperature control.

Technical Data:

Temperature: 64,50C +/- 0,10C
Power supply: 230 Volt AC, 24 Volt DC, 12 VDC

Order Details:

Number Drilling/ Strips	Analog:	Order No.	Number Drilling/ Strips	Digital:	Order No.
max. 4 Strips	12 Volt	4S-12A	10 Drilling	12 Volt	C10B-12D
max. 4 Strips	24 Volt	4S-24A	10 Drilling	24 Volt	C10B-24D
max. 4 Strips	230 Volt	4S-230A	10 Drilling	110 Volt	C10B-110D
10 Drilling	12 Volt	C10B-12A	10 Drilling	230 Volt	C10B-230D
10 Drilling	24 Volt	C10B-24A			
10 Drilling	230 Volt	C10B-230A			

Accessories

Part:	Order No.
Digital reference Thermometer	DK150

