


Stirring Staff PP/ Stirring Cane PVC-white

for Measuring Glass



Application:

for stirring fluids

Technical Details:

Material: Stirring Staff: PP Stirring Cane: PVC white

Order Details:

Part:

Stirring Staff 100/5
Stirring Staff 200/10
Stirring Staff 300/10
Stirring Staff 500/15

Order-No.

RS 100/5
RS 200/10
RS 300/10
RS 500/15

Stirring Cane RST 410/32
Stirring Cane RST 460/46
Stirring Cane RST 540/56
Stirring Cane RST 600/76

RST 410/32
RST 460/46
RST 540/56
RST 600/76

