

Maturation time: 1 min. to 100 h  
Maturation temperature: 0-60°C

Keyboard that optional can be ordered with a pH-control

Pasteurisation time: 1 min. to 100 h  
Pasteurisation temperature: 0°C + 99°C

The housing is built by white polypropylene

The Aerne Analytic Incubators are used for raising bacteria or yeast cultures, e. g. acidification cultures, in the dairying. Optimal results can be reached by programming individual processes. Because of the very good heating isolation of polypropylene the Incubator can be built in a single-walled type

of construction and can anyway be touched by the operator. The Aerne Analytic Incubators are perfect for laboratory uses because of the opportunity of sterile cleaning of the surfaces. All types of Incubators can also be ordered in customized sizes and/ or as tabletop units.