



The Aerne Analytic Magnetic Stirrer Waterbath is basically used for cooking or dissolving of mediums in the microbiology. The instrument is equipped with two magnetic stirrers that can be adjusted continuously. The lid is intended for two bottles with ISO threads for each two or rather five liters, depending on the type. The

waterbath is used with a microprocessor controlling by a keypad. The Magnetic Stirrer Waterbath guarantees an extra high security with its automatically switch off by overheating, sensor break or water shortage. The waterbath is available in series in two different sizes. It can be produced customised in special sizes.