

The secret of perfect heating insulation and space saving design is our mashbath technology made of polypropylene



Standard mashbaths are available with four, eight or twelve mash cups

The Aerne Analytic Mashbaths are used for the implementation of the Kongress- and Hartong method as well as for a variety of free adjustable customized tests. The bath is equipped with five preinstalled programs (out of a possibility of 45), a digital display field, stainless steel containers, exchangeable whisks, a cold water connection with drain cock and draw off taps. The device assures because its made of

polypropylene which makes an ideal heating insulation and furthermore it creates a space saving design. The highest stability of the welds endure variable temperatures. Standard mashbaths are available with four, eight or twelve mash cups. It can also be produced by Aerne Analytic in special sizes. The mashbaths can be used in laboratories because of its possibility to clean the surface sterile.