

Custom insets for various application areas on request



Clock timer with digital timepiece

The Aerne Analytic shaker devices for mojonnier tubes are used to determinate the fat in the milk. The shakers are designed to shake the mojonnier tubes equally. The user defines the time of shaking. The device must be used in a dry laboratory area. The housing is made of PVC and is equipped with a digital

display. An induction motor drives the shaker. Besides the mojonnier shaker Aerne Analytic has a lot of other different products for an effective working with mojonnier tubes in its product range, e. g. centrifuges and/ or heating blocks.