



The Suspended Thermostat Eurostat II of Aerne Analytic is used for heating-up and tempering different fluids in bath tanks. The thermostat is fixed on the rim of vessel with a clamp and is fused against dry heating or overheating by a float switch. The Eurostat II is easy to use by a membrane keyboard. An

automatic shut-off occurs if there is water shortage or a sensor breaks. Furthermore a circulation is provided by a rotor in the medium. The Eurostat II is built by white polypropylene because of its very high resistance to heat. The heating element is of chrome nickel steel.