

Mobile Laboratory Tables

**Application:**

Mobile laboratory benches with continuous safety handle and a three-sided standing board on both working surfaces.

Technical data:

Material working surface:	Polypropylene white, 15 mm material thickness Heat resistance 120°C
Frame:	Aluminium anodised 25×25 mm guide roller with rubber lining (two with locking brake)
Lifting force:	Total 100 kg

Order details:

Dimensions/mm (width x depth)	Measurement/mm Total height	Order No.
700 x 500	750	FBT 70/50/75
700 x 500	900	FBT 70/50/90
800 x 600	750	FBT 80/60/75
800 x 600	900	FBT 80/60/90
1000 x 600	750	FBT 100/60/75
1000 x 600	900	FBT 100/60/90
1200 x 700	750	FBT 120/70/75
1200 x 700	900	FBT 120/70/90
