

# Shaking Shells, Pour-out Caps, Shaking hood

---

## Pour-out Caps for Butyrometers

PVC white

The pour-out cap is placed over the uncorked butyrometers and the rack which are held in place by the cap when poured or rinsed.




### Order Details:

Part	Order No.
12 Test portions	AHB 12
24 Test portions	AHB 24
36 Test portions	AHB 36

## Shaking hood for Butyrometers

PVC white

The shaking hood is placed over supports with sealed butyrometers. The side facing your body is protected by a dropping wall.



### Order Details:

Part	Order No.
12 Test portions	SHB 12
24 Test portions	SHB 24
36 Test portions	SHB 36

## Shaking shells for Butyrometers

Polypropylene (PP) white

The individual butyrometers are placed in the shell. Up to four butyrometers can be shaken, when the cork is held in place simultaneously.



### Order Details:

Part	Order No.
1 Butyrometers	SHÜB 1
2 Butyrometers	SHÜB 2
3 Butyrometers	SHÜB 3
4 Butyrometers	SHÜB 4

---