

Incubator - BRT - Test



Application:

The procedure of inhibitor tests for milk analysis is one of the essential process steps for the quality assurance in the dairying. Test laboratories and farmers use the efficient quick test for their analysis to evaluate the today's high safety standards in this industry. The operating mode of this heating instrument embraces the constant tempering and incubation of the milk test. The Aerne Analytic incubators are equipped with aluminium heating units and an analogue/ digital temperature control.

Order Details: Pre-setted temperature:: 64,5°

Quantity of Slots/Stripes L x B x T 75x9x10 mm	Quantity of Drillings 10,9 mm Drilling	Quantity of Plate(s) 96 Samples	Timer	Volt	Order No. A=Analogue // D=Digital Display
	5			230	5B-230A
	5		Yes	230	5B-230D
	5		Yes	110	5B-110D
	5		Yes	12	5B-12D
	5		Yes	24	5B-24D
	10			230	10B-230A
	10		Yes	230	10B-230D
	10		Yes	110	10B-110D
	10		Yes	12	10B-12D
	10		Yes	24	10B-24D
	24			230	24B-230A
	42			230	42B-230A

	50		Yes	230	50B-230D
	100		Yes	230	100B-230D
2				230	2N00B-230A
2			Yes	230	2N00B-230D
2			Yes	110	2N00B-110D
2			Yes	12	2N00B-12D
2			Yes	24	2N00B-24D
5				230	5N00B-230A
5			Yes	230	5N00B-230D
5			Yes	110	5N00B-110D
5			Yes	12	5N00B-12D
5			Yes	24	5N00B-24D
10			Yes	230	10N00B-230D
10	20		Yes	230	10N20B-230D
20			Yes	230	20N00B-230D
20	20		Yes	230	20N20B-230D
20	40		Yes	230	20N40B-230D
		1	Yes	230	1P00B-230D
		2	Yes	230	2P00B-230D
	20	1	Yes	230	1P20B-230D
	10	1	Yes	230	1P10N-230D