

Incubator CHARM MRL-Test



Application:

The procedure of inhibitor tests for milk analysis is one of the essential process steps for the quality assurance in the dairying. Test laboratories and farmers use the efficient quick test for their analysis to evaluate the today's high safety standards in this industry. The operating mode of this heating instrument embraces the constant tempering and incubation of the milk test. The Aerne Analytic incubators are equipped with aluminium heating units and an analogue/ digital temperature control.

Technical Data:

Control:	Microprocessor
Number Test:	2
Temperature:	55,50C +/- 0,10C
Timer:	8 Minutes with sound signal
Power supply:	230 Volt AC, 24 Volt DC, 12 VDC

Order Details:			
Digital:	Order-No..	Analog:	Order-No..
230 Volt	CT2-230D	230 Volt	CT2-230A
12 Volt	CT2-12D	230 Volt	CT2-230AL
24 Volt	CT2-24D	12 Volt	CT2-12A
110 Volt	CT2-110D	24 Volt	CT2-24A

Accessories

Part:	Order-No..
Transferpette 300 ul Fix	TF300
Pipette tips blau 1000 Stk.	0862
Digital reference Thermometer	DK150