

Incubator Penzym for Chr. Hansen-Test



Application:

The procedure of inhibitor tests for milk analysis is one of the essential process steps for the quality assurance in the dairying. Test laboratories and farmers use the efficient quick test for their analysis to evaluate the today's high safety standards in this industry. The operating mode of this heating instrument embraces the constant tempering and incubation of the milk test. The Aerne Analytic incubators are equipped with aluminium heating units and an analogue/ digital temperature control.

Technical Data:

Control:	Microprozessor
Number Drillings:	10 conic
Temperature:	47,50C +/- 0,10C
Time:	Timer 1: 5 Minutes Timer 2: 10 Minutes
Power Supply:	230 VAC, 24 VDC, 12 VDC
Dimensions (base):	160x100x60 mm
Weight:	0,8 kg

Order Details:			
Digital:	Order-No.	Analog:	Order-No.
230 Volt	P10-230D	230 Volt	P10-230A
12 Volt	P10-12D	12 Volt	P10-12A
24 Volt	P10-24D	24 Volt	P10-24A

Accessories

Part:	Order-No.
-------	-----------

Transferpette 10 ul Fix	TF10
Transferpette 10 ul Fix	TF50
Transferpette 100 ul Fix	TF100
Reactors 500 Pcs.	0760
Pipette tips yellow 1000 Pcs.	0132
Pipette tips white 200 Pcs.	1350
Digital reference Thermometer	DK150
Stand for Reactors	ST20