


Packaging Tester (Leakage Tester)

VT250 / VT350 / VT-Q



Application:

The Packaging Tester is established to test nearly all packagings under vacuum. The Aerne Analytic Vacuum Tester is developed with a digital pressure gauge and an automatic ventilation. This random analysis is exactly what the today's quality standards in the packaging industry require. This user-friendly device is perfect for laboratory use because the surfaces can be cleaned sterile. The pressure stability assures with its special drafted and well-tried components. To use with different packaging sizes there are a multitude of Vacuum testers in the Aerne Analytic range of products.

Custom-made Leakage Testers that have customised specifications are possible to produce by Aerne Analytic as well.

Order Details:

Part:

VT 250 (analog)

- Fittings: Vakuummeter with Precision regulator 0 to 1000 mbar, Water supply connection, bleeder screw
- Immersion plate revolving and height-adjustable. Self-closing lid. To connect diaphragm pumps (not included in the scope of supply).
- Material test Chamber: Plexi glass clear
- Material Housing: PVC white

VT 350 (analog)

- Fittings: Vakuummeter with Precision regulator 0 to 1000 mbar, Water supply connection, bleeder screw
- Immersion plate revolving and height-adjustable. Self-closing lid. To connect diaphragm pumps (not included in the scope of supply).
- Material test Chamber: Plexi glass clear
- Material Housing: PVC white

Diameter	Height	Order-No.
200 mm	250 mm	VT 250
300 mm	250 mm	VT 350

VT 250D (digital)

- Jumo Controller for Vacuum: 0-500 mbar
- digital Display
- Incl. pre-defined Vacuum Ramp (to self conform)
- round Container
- Material test Chamber: Plexi glass clear
- Material Housing: PVC white
- Power Supply : 230V 50/60 Hz

200 mm 250 mm VT 250D

VT 350D (digital)

- Jumo Controller for Vacuum: 0-500 mbar
- digital Display
- Incl. pre-defined Vacuum Ramp (to self conform)
- round Container
- Material test Chamber: Plexi glass clear
- Material Housing: PVC white
- Power Supply : 230V 50/60 Hz

300 mm 250 mm VT 350D

Dimensions
VT QD

- Jumo Controller for Vacuum: 0-500 mbar
- digital Display
- Incl. pre-defined Vacuum Ramp (to self conform)
- square Container
- Material test Chamber: Plexi glass clear
- Material Housing: PVC white
- Power Supply : 230V 50/60 Hz
- **Customized special sizes on request!**

540 x 440 x 320 mm VT QD

Flyer:

