

Autoclave



Application:

Autoclave for sterilization up to 95°C

Technical Data:

Casing:	Material Polypropylene white enclosed lid
Temperature controller:	0 + 100°C +/- 0,5°C
Heating:	Chromium nickel steel-Tube heating element
Safety:	Switches off when water is lacking. Overheating protection.
Power supply:	230 VAC

Order Details Model 1 (Payload Dimensions/mm: LxBxH 400x400x450):

Part

Model with Timer 120 Min. and Temperature-Controller 0+100°C
Aluminium Support, Square tube 25 mm

Order No.

DTM 1
AUB

**Order Details Model 2 (Payload Dimensions/mm: LxBxH
300x300x520):****Part**

Model with Timer 120 Min. and Temperature-Controller 0+100°C
Aluminium Support, Square tube 25 mm

Order No.

DTM 2
AUB

Insert basket of stainless steel:**Part :**

Basket

Order No.

EK
