


Revolving Bath External Type 1P/ Type 2P

**Application:**

The Revolving Baths are utilised for the heating of external containers or baths in laboratories.

Technical Data:

Casing:	Polypropylene (PP) white Material strength 10 mm
Heating element:	Chromium nickel steel 1000 watt / 2000 watt variable 10°C to 99°C +/- 0.3°C (depending on model)
Control:	minimum temperature depending on ambient temperature
Setting range:	10°C bis 99°C
Pump performance Type 1P:	12 l/min. with inflow and outflow tube
Pump performance Type 2P:	2 x12 l/min. individually Connectable with 2 inflow and outflow tubes each
connections:	1 inflow and outflow tube each
Safety:	Overheating protection and level protection
Power Supply:	230 VAC

Order Details:

Type	Payload dimensions cm LxWxH	Numbers of Pumps	Order-No.
P1	24x20x17	1	UWT 1
P2	29x25x17	2	UWT 2
